

## Shower Stool with back



**Brand:** Aidapt  
**Product Code:** VB540S  
**Availability:** 25  
**Weight:** 2.90kg

**Price: £63.00**

### Short Description

The Shower Stool with Back is adjustable in height to suit most users and has integral handles in the seat and back to aid support and transfer, with non-slip feet for safe use in wet conditions.

### Description

A high quality shower stool with non-slip feet designed for everyday use. The Shower Stool with Back is adjustable in height to suit most users, with clearly-numbered height settings, and has integral handles in the seat and back to aid support and transfer. The lightweight frame means it can be removed with ease enabling multiple use of the shower area.

- Perforated top and back for easy drainage
- 8 numbered height adjustment settings
- Tubular steel frame featuring integral grips for extra support
- Non-slip feet for safe use in wet conditions
- Back Rest (mm): 390x240
- Colour: WHITE
- Depth (mm): 420
- Footprint (WxD) Min/Max (mm): 490x430
- Height (mm): 700-928
- Maximum User Weight (kg): 120

- Maximum User Weight (stones): 18.75
- Net weight (kg): 2.9
- Overall Height (mm): 690-865
- Product Dimensions (mm): 700-928x500x420
- Seat Depth (mm): 290
- Seat Height (mm): 390-565
- Seat Width (mm): 490
- Width (mm): 500

**VAT RELIEF APPLIED**  
[Read more about VAT Exemption...](#)

## Product Gallery

