

Compact folding shower stool



Brand: Aidapt
Product Code: VB541
Availability: 25
Weight: 3.30kg

Price: £49.99

Short Description

This space-saving Folding Shower stool/chair is an ideal solution if space is at a premium and folds easily for storage.

Description

This space-saving Folding Shower Chair is an ideal solution if space is at a premium. The steel tube frame is powder coated to prevent corrosion and the anti-slip rubber feet aid safety and stability. The blow moulded HDPE seat has been designed with built-in handles and drainage holes. The chair folds easily for storage.

Specifications

- Space saving folding shower chair allowing independent bathing
- Maximum User Weight Limit 100kg (16 stone)
- Anti-Slip Rubber Feet for Stability
- Powder-coated Steel for durability
- Moulded HDPE Seat
- Colour: White
- Depth (mm): 300
- Floor to Arm (mm): 505
- Footprint (WxD) Min/Max (mm): 430x390

- Height (mm): 470
- Height of Handle from Floor (mm): 505
- Height To Top Of Back From Floor: 470
- Maximum User Weight (kg): 100
- Maximum User Weight (stones): 16
- Net weight (kg): 3.3
- Overall Depth (mm): 300
- Overall Height (mm): 470
- Product Dimensions (mm): 500x300x470
- Seat Depth (mm): 290
- Seat Height (mm): 470
- Seat Width (mm): 490
- Width (mm): 500

[Read more about VAT Exemption](#)

Product Gallery

