

## Round adjustable shower stool



**Brand:** Aidapt  
**Product Code:** VB511A  
**Availability:** 25  
**Weight:** 1.00kg

**Price: £34.65**

### **Short Description**

Multi purpose, versatile, adjustable aluminium stool.

### **Description**

This versatile, multi-purpose aluminium stool can be used in a whole host of different settings around the home, in the garden or playroom and even when camping. With a comfortable moulded plastic seat top, this affordable and contemporary design is durable. Constructed from a robust aluminium frame, it's lightweight, waterproof, rust-resistant and height adjustable with non-slip ferrules for safety and comfort. Small in size when the legs are compacted and stackable - for easy storage and transporting. Supplied flat packed, with simple home assembly required, can be used indoors or outdoors and suitable for adults or children.

### **Specifications**

- Multi-purpose stool, can be used both indoors and outdoors
- Maximum user weight limit 120kgs (18.9 stones)
- Durable, hard-wearing materials with rust-proof aluminium frame
- Increased grip ensures stability on damp or wet floors
- Height adjustable legs and non-slip ferrules
- Stackable for space-saving storage and transportation

- Weight (kg) 1
- Seat Height(mm) 350 x 525
- Seat Depth (mm) 315
- Seat Width 315
- Foot print (wxd) (mm) 330 x 330
- Maximum user weight 120 kg (19 stone)

[Read more about VAT Exemption...](#)