

## Pembury Height Adjustable Corner Shower Stool



**Brand:** Aidapt  
**Product Code:** VB636  
**Availability:** 25  
**Weight:** 2.00kg

**Price: £66.92**

### Short Description

Featuring fine leg adjustment to maximise stability and safety, and a clip-on/clip-off seat for easy cleaning.

### Description

The Pembury Height Adjustable Corner Shower Stool is an excellent solution when space is at a premium. It features fine leg adjustment to maximise stability and safety when in use and a clip-on/clip-off seat for easy cleaning. This compact, corner shower stool is manufactured from rust free aluminium for improved longevity and service life..

### Specifications

- Corner shower stool is ideal where space is at a premium
- Robust clip-on/ clip-off seat for easy drainage
- Simple height adjustment leg to suit most users
- Manufactured from rust-free aluminium
- 159 kg (25 stone) weight limit
- Colour: White
- Depth (mm): 380

- Footprint (WxD) Min/Max (mm): 350-550 x 380
- Height (mm): 470-620
- Maximum User Weight (kg): 159
- Maximum User Weight (stones): 25
- Net weight (kg): 2
- Seat Depth (mm): 265
- Seat Height (mm): 470-620
- Seat Width (mm): 410

 [Read more about VAT Exemption...](#)