

Coniston Rotating Bather



Brand: Aidapt
Product Code: VR130
Availability: 15
Weight: 6.10kg

Price: £181.50

Short Description

The Coniston Rotating Bather, high quality bath seat, is width adjustable and the seat can turn 90 degrees in either direction.

Description

The Coniston Rotating Bather is a high quality rotating bath seat. It is width adjustable, fitting not only standard baths, but also extra wide baths (up to 27/28"). The Coniston Bath Seat has a pivoting base that enables the seat to turn 90 degrees in either direction, it can be locked in position for maximum stability. The Coniston Bather offers maximum stability and comfort to the user making it an ideal alternative to many bath boards, it also has the benefit of width adjustable arms. No special fittings are needed to install and the Coniston Rotating Bather can be removed in seconds.

Specifications

- Four positive seat locking positions
- Width adjustable bars, suitable for 26, 27 and 28 wide baths
- Moulded epoxy coated plastic one-piece seat which conforms to MHRA guidelines
- Easy to clean, all alloy frame for longevity suitable for industrial cleaners
- 150 kg (24 stone) weight limit

- Colour: White
- Configuration: Aluminum
- Maximum User Weight (kg): 150
- Maximum User Weight (stones): 24
- Net weight (kg): 6.1
- Overall Depth (mm): 610
- Product Dimensions (mm): 610x695x610
- Seat Height (mm): 60
- Seat Width (mm): 450
- Weight Limit (kg): 150
- Width Adjustment (mm): 695-750
- Width between Arms (mm): 570-670

SAFETY NOTICE

It is recommended that a minimum of 30mm (1¼) of each arm must rest securely on each side of the bath and the rim of the bath be level.

[Read more about VAT Exemption...](#)

