

Medina plastic bath seat



Brand: Aidapt
Product Code: VR127
Availability: 20
Weight: 2.20kg

Price: £52.80

Short Description

The Medina Bath Seat provides extra height enabling easy access in and out of the bath.

Description

The Medina Bath Seat provides extra height enabling easy access in and out of the bath. The Medina Bath Seat is constructed from high quality, fully reinforced slatted plastic and meets MHRA safety guidelines. The bath seat is secured to the bath with 4 large suction caps which can be removed by pulling on the integral tabs. Available in 3 heights to suit all users.

Specifications

- Modern easy to clean surface
- Slatted for easy drainage, conforms to MHRA guidelines
- Finished in soft grey
- 190 kg (30 stone) weight limit
- Retail boxed

Code

Colour

Depth (mm)

VR127A

White and grey

290

VR128A

White and grey

290

Height (mm)	152	305
Maximum User Weight (kg)	190	190
Maximum User Weight (stones)	30	30
Net weight (kg)	2.2	
Overall Depth (mm)	290	290
Packing	Retail Box	Retail Box
Product Dimensions (mm)	152x450x290	305x450x290
Seat Height (mm)	152	305
Seat Width (mm)	450	450
Size	152 mm (6")	304 mm (12")
Weight Limit (kg)	190	190
Width (mm)	450	450

[Read more about VAT Exemption](#)

Product Gallery

