

## Pride Apex Aluminum



**Brand:** Pride  
**Product Code:** ALUMALITE  
**Availability:** 30  
**Weight:** 28.80kg

**Price: £2,299.00**

### Short Description

The Pride Apex Aluminum is the lightest aluminum travel scooter on the market.

### Description

The Pride Aluminium scooter is the lightest aluminum travel scooter on the market! Experience the convenience of a lightweight, compact travel scooter. Easy to put into the boot, with its heaviest piece weighing only 11 kg. The Apex Alomalite, with its aluminum frame, weighs only 28.8kg in total making it our lightest four-wheeled scooter yet! **Specifications:**

- Weight Capacity: 18 stone 13 lbs (120 kg)
- Maximum Speed: 4 mph (6Kph)
- Ground Clearance: 32 mm (1.25")
- Turning Radius: 1385 mm (54.5")
- Overall Length: 1022 mm (40.25")
- Overall Width: 490 mm (19.25")
- Drive Tyres: 203 mm (8") front & rear
- Total Weight Without Batteries: 28.8 kg (63.4 lbs)
- Heaviest Piece When Disassembled: 11.2 kg - main frame
- Dual Breaking System Regenerative & electromechanical
- Battery Requirements: 12 Ah Lithium Battery
- Battery Weight: 3.5 kgs (8 lbs)

- Battery Charger: 2 amp off-board

**Colour: Red, Blue, Black**

[Read more about VAT Exemption.](#)

## Product Gallery

