Pride Apex Rapid



Brand: Pride

Product Code: APEXRAP

Availability: 60 Weight: 39.60kg

Price: £1,789.00

Short Description

Pride Apex Rapid is featuring all round suspension mobility scooter.

Description

Pride Apex Rapid? with luxury and convenient features like a wraparound delta tiller with a high level charging port, LED kerb light and a storage basket together with a compact, easy to transport design.

Sleek styling complements the Apex Rapid's abundant features making it an exceptional travel mobility choice.

Specifications:

Battery Weight (each) Weight: 18 AH: 13 kg (29 lbs.)

Weight Capacity 21 Stone 6 lbs (136 kg)

Maximum Speed 4 mph (6 km/h)

Ground Clearance 3.8 cm (1.5")

Overall Length 104.1 cm (41')

Overall Length 104.1 cm (41')
Overall Width 52.1 cm (20.5")

Turning Radius 116.2 cm (50.5")

Standard Controls LED Display

Heaviest Piece (when disassembled) 16.6 kg (36.6 lbs.) front section

Total Weight (without batteries)

Suspension

Battery Requirements

Seating

Tyres

39.6 kg (87.4 lbs)

Front and rear independent CTS

Suspension

Two 12-volt, deep cycle (standard) Size:

2 x 18 AH

Foldable, compact, vinyl-covered

molded plastic

Front: 5 cm x 17.8 cm (2" x 7") solid

Rear: 6.35 cm x 20.3 cm (2.5" x 8")

solid

Read more about VAT Exemption...

Product Gallery









