

Sock and Stocking Aid



Brand: Aidapt
Product Code: VM944
Availability: 30
Weight: 0.00kg

Price: £9.90

Short Description

Sock and Stocking Aid helps people to put on socks or stockings without bending down to the feet.

Description

The Sock and Stocking Dressing Aid makes putting on socks and stockings easy. Large loop handles at the end of strong 32 inch straps make this a popular dressing aid for arthritis sufferers. The nylon lining of the sock and stocking aid helps the foot to slide easily through the garment.

- Handy aid that helps the user put on socks reducing the need to bend
- Can use with majority of socks and stockings

- Depth (mm): 180
- Height (mm): Soft
- Width (mm): 230
- Product Dimensions (mm): 230x180x10
- Net weight (kg): 0.1
- Colour: White and Blue

Product Gallery

